



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Skin Closure Strips Specifications

Catalogue No. /Reference No. /Products Code	JC05022
Material	Porous non woven fabric; Pressure-sensitive hypoallergenic adhesive; Reinforced with polyester filaments for added strength; Backing paper
Color	White Color
Size	3x75mm, 6x38mm, 6x75mm, 6x100mm, 12x100mm, etc
Package	1 or 2 Pcs/Sterile Pouch
Shelf Life	5 Years
Sterile or Not	Sterile
Rules & Classification	Rule 4, Class I Sterile





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Skin Closure Strips Specifications

Reference No/ Product Code	Description & Size	Size	Unit Package	Outside/Car- ton Package	Carton Size
JC050221	Wound Closure Skip Strip, 3x75mm, 5 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	3x75mm	5 strips/pouch, 50 pouches/box	20 boxes/carton	30.5x15x14CM
JC050222	Wound Closure Skip Strip, 6x38mm, 6 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	6x38mm	6 strips/pouch, 50 pouches/box	20 boxes/carton	30.5x15x14CM
JC050223	Wound Closure Skip Strip, 6x75mm, 3 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	6x75mm	3 strips/pouch, 50 pouches/box	20 boxes/carton	30.5x15x14CM
JC050224	Wound Closure Skip Strip, 6x100mm, 5 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	6x100mm	5 strips/pouch, 50 pouches/box	20 boxes/carton	36x32x18CM
JC050225	Wound Closure Skip Strip, 6x100mm, 10 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	6x100mm	10 strips/pouch, 50 pouches/box	20 boxes/carton	36x32x18CM
JC050226	Wound Closure Skip Strip, 12x100mm, 6 Strips/sterile pouch, 50pouches/box, 20 Boxes/Carton, EO Sterilization.	12x100mm	6 strips/pouch, 50 pouches/box	20 boxes/carton	36x32x18CM